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
EN

IMPORTANT! READ CAREFULLY AND KEEP FOR FUTURE REFERENCE.

I - WARNING

- Please read the instructions carefully before using the digital thermometer.
- Risk of choking. The thermometer cap and battery can be fatal if swallowed. Do not allow children to use this device without adult supervision.
- Do not use this thermometer in the ear. It is designed to take the temperature in the mouth, rectum or under the armpit.
- Do not expose the battery to a high heat source as it may explode.
- Note: Use of the cover on the thermal probe may result in a deviation of 0.1°C (0.2°F) from the actual temperature.
- When not using the thermometer for an extended period of time, remove the battery.
- Using temperature readings for self-diagnosis is dangerous. Consult your physician for interpretation of your readings. Self-diagnosis can lead to worsening of disease symptoms.
- Do not attempt to take the temperature when the thermometer is wet. This may result in inaccurate readings.
- Do not bite the thermometer. This could break the thermometer and/or cause injury.
- Do not attempt to disassemble or repair the thermometer. This may result in inaccurate readings.
- After each use, disinfect the thermometer, especially if it is to be used by more than one person.
- Never force the thermometer into the rectum. In case of pain, remove the thermometer and discontinue taking the temperature. Follow this instruction to avoid risk of injury.
- Do not put the thermometer in the mouth when it has been used to take a rectal temperature.
- Do not use the thermometer to take the oral temperature of a child two years of age or younger.
- When the thermometer has been exposed/stored at a temperature below 5°C (41°F) and above 40°C (104°F), the temperature should be taken at a temperature between 5°C (41°F) and 40°C (104°F). ~ 104°F), leave it at room temperature for about 15 minutes before using it.

• PLEASE READ THESE INSTRUCTIONS CAREFULLY BEFORE USING THE THERMOMETER:


- This digital thermometer provides a quick and extremely accurate reading of body temperature. It is designed to measure the temperature body in the mouth, rectum or under an armpit on a daily basis. It can be reused in the clinic or at home with people of all ages. To better understand how it works, and to benefit from years of reliable use, please read all instructions before using it.
- This thermometer complies with the following standards : EN 12470-3 Clinical thermometers - Part 3: Performance of compact electrical thermometers (with or without prediction) with maximum device. ISO 80601-2-56 Medical electrical equipment - Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for the measurement of body temperature. EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance.
– Collateral Standard: Requirements for medical electrical equipment and electromedical systems used in the home and in accordance with requirements safety standards EN 60601-1-2 (EMC), IEC/EN 60601-1. The manufacturer is ISO 13485 certified.

II - CONTENT

- 1 Thermometer. 1 User manual. 1 Storage case.

III - PRODUCT IMAGE

Figure 1

IV - WARNING

*The performance of the device may degrade in the following cases :

- Use outside the thermal and hygrometric values stipulated by the manufacturer.
- Storage outside the thermal and hygrometric values stipulated by the manufacturer.
- Mechanical shock (e.g. drop test) or damaged probe.
- Measured body temperature below ambient temperature.

*Portable and mobile RF communications may interfere with the operation of the device. This device requires special precautions to protect it from electromagnetic interference (EMI) in accordance with the EMI information provided in the manufacturer's documentation.

*Keep the thermometer away from magnetic resonance equipment.

V - EXPLANATION OF SYMBOLS

	Applied part Type BF		Manufacturer
	Consult the documents provided by the manufacturer	55°C -20°C 	Storage and transport Thermal limits : -20°C ~ 55°C (-4°F ~ 131°F)
	Authorized European Representative	LOT :	Lot Code

EN

VI - FEATURES

Type :	Digital thermometer (without prediction)
Measuring range:	32.0°C ~ 42.9°C (90.0°F - 109.9°F) (°C /°F selected by the manufacturer)
Accuracy :	±0,1°C (±0,2°F) à 35,5°C ~ 42,0°C (95,9°F ~ 107,6°F) à 18°C ~ 28°C (64,4°F ~ 82,4°F) Ambient temperature range ±0.2°C (±0.4°F) for other measurements and operating temperature range
Operational mode:	Direct mode
Display :	3.5-digit LCD display
Memory :	To save the last measured value
Stack :	One 1.5 V DC button cell battery (LR41 or SR41, UCC 392)
Autonomy :	Approximately 200 hours of continuous operation or one year with three measurements per day.
Dimensions :	13,1cm x 2,0cm x 1,2cm (L x larg. x H)
Weight :	Approx. 10 grams with battery
Expected lifetime:	Three years
Ambient operating temperature range:	Temperature: 5°C ~ 40°C (41°F ~ 104°F) Relative humidity: 15% ~ 95%RH Atmospheric pressure: 700hPa ~ 1060hPa
Storage and transportation	Condition: Temperature : -20°C ~ 55°C (-4°F ~ 131°F) Relative humidity: 15% ~ 95%RH Atmospheric pressure: 700hPa ~ 1060hPa
Index of protection:	IP 27
Classification :	Type BF

°C/°F SWITCH

Temperatures can be measured in degrees Celsius or Fahrenheit (°C/°F; displayed in the upper right corner of the LCD screen). When the thermometer is off, press and hold the On/Off button for approximately 2 seconds to change the unit of measure.

VII - INSTRUCTIONS FOR USE


Figure 2

5. Press the On/Off button located next to the LCD screen. A tone sounds and the display shows **188.8%** followed by the last measured temperature. After the self-checking procedure, the thermometer is in temperature-taking mode. Position the thermometer in the

place chosen for temperature taking (mouth, rectum, armpit).

- a) **Taking oral temperature:** Place the thermometer under the tongue in the position shown in Figure 2. Close your mouth and breathe in normally and regularly through your nose to avoid disturbance of the measurement by exhaled/inhaled air. A normal body temperature is between 35.7°C and 37.3°C (96.3°F and 99.1°F).
- b) **Rectal temperature taking:** Lubricate the silver end of the catheter with Vaseline to facilitate insertion. Gently insert the thermometer probe into the rectum for a length of one centimeter (less than 0.5 inch). A normal body temperature is between 36.2°C and 37.7°C (97.2°F and 99.9°F).
- c) **Taking the temperature of an armpit:** Wipe the armpit dry. Place the probe in the armpit and squeeze the arm firmly on the side of the body. From a medical point of view, this method always gives inaccurate readings. Therefore, it should not be used if specific temperatures are required. A normal body temperature is between 35.2°C and 36.7°C (95.4°F and 98.1°F).

6. The Degree symbol flashes during the entire temperature measurement. When it stops flashing, an alarm sounds for about 10 seconds. Simultaneously, the measured temperature is displayed on the LCD. The minimum duration of the temperature measurement up to the alarm (beep) must be observed without exception. The measurement continues even after the alarm signaling the end of the measurement. To obtain a better measurement of the body temperature, we recommend leaving the probe in the mouth or rectum for at least 2 minutes, and for 5 minutes in the armpit, regardless of the alarm. A minimum measuring interval of 30 seconds must be observed.

7. To extend the battery life, press the ON/OFF button to turn the thermometer off after the temperature reading is completed. If you do not press this button, the thermometer will automatically turn off after ten minutes.

*Note : Normally, the alarm sounds "Bi-Bi-Bi-Bi-" and accelerates when the temperature reaches 37,8°C Or more. The alarm then sounds "Bi-Bi-Bi----- Bi-Bi-Bi--Bi-Bi-Bi".

VIII - TROUBLESHOOTING

Error Message	Problem	Solution
Lo	Temperature displayed is below 32.0°C (90.0°F)	Turn off the thermometer. Wait one minute and measure the temperature again with direct contact and after a sufficiently long interval.
Hi	The temperature displayed is higher at 42.9°C (109.9°F)	Turn off the thermometer. Wait one minute and measure the temperature again with direct contact and after a sufficiently long interval.
Err	The system is not working properly.	Remove the battery. Wait one minute, then turn the thermometer back on. If the error message persists, contact the supplier.
	Battery empty : The battery icon flashes but the capacity is not measurable.	Replace the battery.

IX - BATTERY REPLACEMENT

- Replace the battery when  is display in the corner bottom right of the LCD screen.
- Remove the battery cover as shown in Figure 3.
- Gently pull the plastic circuit board out of the plastic circuit board by about 1 cm with the battery compartment (just under 0.5 inch, see Figure 4).
- Use a sharp object, such as a ballpoint pen, to remove the empty battery. Dispose of the battery by respecting the standards and the law. Install a new 1.5V DC button cell battery type LR41 or SR41, UCC392 or equivalent. Check the polarities: the + sign must face you. (See Figure 5).
- Reinsert the battery compartment in position and replace the battery cover.


Figure 3

Figure 4

Figure 5

X - CLEANING AND DISINFECTION

Wipe the thermometer with a clean, soft cloth.

For tougher stains, wipe the thermometer with a cloth dampened with water or a solution of a mild detergent of neutral detergent. Rinse it thoroughly. Finally, wipe with a soft, dry cloth.

To disinfect it, use an agent based on isopropyl alcohol or 75% ethanol.

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- Do not attempt to disinfect the thermometer probe (tip) by immersing it in alcohol or hot water (above 50°C/ 122°F).
- Do not use an ultrasonic system to clean the thermometer.

XI - CALIBRATION

The thermometer is calibrated by its manufacturer. If the thermometer is used in accordance with the instructions in this manual, periodic adjustment is not necessary.

Nevertheless, we recommend checking the calibration every two years, or whenever the clinical accuracy of the thermometer is in question. Turn on the thermometer and insert the probe into a water bath. Then check the accuracy of the thermometer. Send the entire device to the supplier or manufacturer.

The above recommendations do not replace legal requirements. The user must always comply with the legal requirements applicable to measurement control, functionality and accuracy of measurements, which are defined by laws, directives or ordinances, concerning the use of the thermometer.

XII - LIMITED WARRANTY

The thermometer has a one-year warranty from the date of purchase. If the thermometer does not work properly due to defective components or a manufacturing defect, we will repair or replace it free of charge. All components are covered by this warranty, excluding the battery. The warranty does not cover damage to the thermometer caused by improper use or handling. To receive warranty service, you must present the original proof of purchase or a copy provided by the original seller.

Disposal of the product and empty batteries must comply with national regulations for the disposal of electronic products.

The product complies with the requirements of the directive MDD 93/42/0197 EEC.

The regulatory body is identified by the number « 0197 ».

EN

XIII - INFORMATION ON ELECTROMAGNETIC COMPATIBILITY(EMI)

This device complies with the EMI requirements of the international standard IEC 60601-1-2. The requirements are met in accordance with the conditions stipulated in the table below. This device is an electrical medical device, which must be used with special precautions to avoid electro-magnetic interference as indicated in its instructions for use. Portable and mobile RF communications equipment may interfere with the operation of this device. The use of this thermometer with non-approved accessories may adversely affect its operation and electromagnetic compatibility. This device should not be used in the direct vicinity of other electrical equipment or between two electrical devices.

Table 1

Information and manufacturer's declaration on electromagnetic emissions		
The device is designed to be used in the following environments electromagnetic specified below. It is the responsibility of the customer or user to verify that the device is always used in such an environment.		
Emissions testing	Compliance	Advice on the electromagnetic environment
RF emissions CISPR 11	Group 1	The device uses only RF energy for its internal operation. Therefore, its RF emissions are very low and are unlikely to cause interference to nearby electronic equipment.
RF emissions CISPR 11	Class B	The device can be used in any building or house, including domestic establishments and those directly connected to a public low-voltage power grid that supplies private residential buildings.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/transient emissions IEC 61000-3-3	N/A	

Table 2

Information and manufacturer's declaration on electromagnetic immunity of the manufacturer			
The device is designed for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to verify that the device is always used in such an environment.			
Immunity test	Test level IEC 60601	Level of compliance	Advice on the electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors must be wood, concrete or tile. If the floor is made of synthetic material, the humidity relative must be at least equal to or greater than 30%.
Electrostatic transients/bursts IEC 61000-4-4	Power supply lines ±2 kV Input/output lines ±1 kV	N/A	
Surtension IEC 61000-4-5	Differential mode ± 1 kV Common mode ± 2 kV	N/A	
Chute de tension, courtes interruptions, variations de tension sur les lignes d'alimentation IEC 61000-4-11	< 5% de l'UT (baisse >95% de l'UT) pour 0,5 cycle 40% de l'UT (baisse de 60% de l'UT) pour 5 cycles 70% de l'UT (baisse de 30% de l'UT) pour 25 cycles <5% de l'UT (baisse >95% de l'UT) pour 5 cycles	N/A	
Electromagnetic frequency field (50/60 Hz) IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz	Frequency fields the magnetic properties of the power supply must be compatible with the standard requirements of a commercial environment or typical hospital.

EN

Table 3

Information and manufacturer's declaration on electromagnetic immunity of the manufacturer			
The device is designed to be used in the following environments electromagnetic specified below. It is the responsibility of the customer or user to verify that the device is always used in such an environment.			
Immunity test	Test level IEC 60601	Level of compliance	Advice on the electromagnetic environment
Pipes RF IEC 61000-4-6	3 Vrms 150kHz à 80MHz	N/A	<p>Portable and mobile RF communications equipment must be kept separate from the thermometer at all points, including cables, and at a distance of recommended separation calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separationn distance</p> $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ <p>800 MHz to 2.7 MHz</p> <p>P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer's data, and d is the recommended separation distance in metres (m). The field strengths of fixed RF transmitters, as determined by an assessment of electromagnetic emissions of the site must be below the compliance level in each frequency range. Interference may occur near an equipment identified by the following symbol :</p> 
Radiated RF IEC 61000-4-3	10 V/m 80MHz à 2.7Ghz	10 V/m	
RF wireless communication equipment IEC 61000- 4-3	380MHz, 27V/m	380MHz, 27V/m	
	450MHz, 28V/m	450MHz, 28V/m	
	710MHz, 745MHz, 780MHz 9V/m	710MHz, 745MHz, 780MHz 9V/m	
	810MHz, 870MHz, 930MHz 28V/m	810MHz, 870MHz, 930MHz 28V/m	
	1720MHz, 1845MHz, 1970MHz 28V/m	1720MHz, 1845MHz, 1970MHz 28V/m	
	2450MHz, 28V/m	2450MHz, 28V/m	
	5240MHz, 5500MHz, 5785MHz 9V/m	5240MHz, 5500MHz, 5785MHz 9V/m	

Table 4

Recommended separation distances between the thermometer and portable and mobile RF communication equipment		
<p>The thermometer is designed to be used in a radiant electromagnetic environment with controlled interference.</p> <p>The customer or the user of the thermometer can prevent electromagnetic interference by maintaining an minimum distance between the thermometer and portable and mobile RF communication equipment (transmitters) by respecting the recommendations below, and according to the maximum output power of the communication equipment.</p>		
Maximum transmitter output power - W	Separation distance according to transmitter frequency - m	
	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_3} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{7}{E_3} \right] \sqrt{P}$
0,01	0,12	0,23
0,1	0,38	0,73
1	1,2	2,3
10	3,8	7,3
100	12	23
<p>For transmitters whose maximum output power is not indicated above, the separation distance (d) in meters (m) can be estimated using the equation applicable to the transmitter frequency, where P is the maximum output power of the transmitter measured in Watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range must be applied.</p> <p>NOTE 2: These recommendations may not apply to all situations. Electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.</p>		

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JOYTECH HEALTHCARE Co. Ltd.
Shanghai International Holding Corp. GmbH (Europe)